

# NEW Steelco Tip Descaler

Effective automated scaler tips processing solution for Steelco washer disinfectors



Steelco provides a solution to achieve professional validated cleaning and disinfection for different brands and sizes of scaler tips.



Tip Descaler allows fast flexible loading and removal of up to three scaler tips. Just replace one or more of the hand piece adaptors in the wash cart of your Steelco washer disinfectors.



Tip Descaler is made of high quality flexible silicone for extended life.



Tip Descaler holder is compatible with all Steelco wash carts for cleaning dental handpieces.



A fine mesh filter reduces the risk of debris clogging the jet of the scaler tip.

Headquarters

**STEELCO S.p.A.**

Via Balegante, 27  
31039 Riese Pio X (TV) - Italy  
Ph. +39 0423 7561  
Fax +39 0423 755528  
e-mail: [steelco@steelcospa.com](mailto:steelco@steelcospa.com)

[www.steelcospa.com](http://www.steelcospa.com)

Branches

**STEELCO HUNGARY**

Temesvár u., 19-21  
1116 Budapest - Hungary  
e-mail: [steelco@t-online.hu](mailto:steelco@t-online.hu)

**STEELCO DENMARK**

Lottenborgvej, 24  
2800 Kgs. Lyngby - Denmark  
[steelco-denmark@steelcospa.com](mailto:steelco-denmark@steelcospa.com)

**STEELCO USA**

999 Stinson Way - Unit 307  
West Palm Beach, FL 33411 USA  
e-mail: [steelco-usa@steelcospa.com](mailto:steelco-usa@steelcospa.com)

**STEELCO INDIA**

Plot No. 877, Sector 40,  
Gurugram - 122001 - Haryana - India  
[steelco-india@steelcospa.com](mailto:steelco-india@steelcospa.com)

**STEELCO ASIA**

6th floor, Suite 16, IOI Business Park - Persiaran Puchong Jaya Selatan  
Bandar Puchong Jaya - 47100 Puchong, Selangor - Malaysia  
e-mail: [steelco-asia@steelcospa.com](mailto:steelco-asia@steelcospa.com)

**STEELCO MEXICO**

Calle Guadalupe, 16 - Colonia Pueblo Nuevo,  
Delegación La Magdalena Contreras, 10640 CDMX - Mexico  
[steelco-mexico@steelcospa.com](mailto:steelco-mexico@steelcospa.com)