

NEW Steelco TTA 18 - TTA 23
Compact Class B sterilizers
for dental and medical clinics

- **Safe & Easy**
- **Flexible loading rack**
- **Designed to fit into your facility**



Volume chamber
 18 lt / 24 lt

Designed for the fast and efficient safe sterilization of solid or hollow wrapped and unwrapped instruments and porous loads.

Steelco's TTA 18 and TTA 23 models are **configured for either steam generation in the chamber or separate steam generator within the autoclave body**. Can be connected to mains water supply or operated as a stand-alone unit with manual water filling. A wide choice of accessories is available.



Pouch loading insert



Aquafilter for the production of deionized water



Aquadist for distilled water



PS 2001 heat sealer pouches

NEW Steelco Tip Descaler
Effective automated scaler tips processing solution for Steelco washer disinfectors

Steelco provides a solution to achieve professional validated cleaning and disinfection for different brands and sizes of scaler tips.



Tip Descaler allows fast flexible loading and removal of up to three scaler tips. Just replace one or more of the hand piece adaptors in the wash cart of your Steelco washer disinfectant.

Check it out!



37th INTERNATIONAL DENTAL SHOW

Join us at IDS in Cologne from March 21st to 25th
 Steelco looks forward to introduce you to these exciting new products. **Hall 10.1 - Stand D058 E059**