

## What is the level of treatment required for your endoscopes?

Now you have the freedom to choose the most suitable and safe process!



### DOUBLE ISO 15883-4 & ISO 14937 COMBINATION COMPLIANCE



Prevention of multidrug-resistant health care associated infections (HAI) has become a priority issue and great challenge to clinicians.

For better control of HAI and in line with customers needs and endoscope reprocessing latest requirements Steelco S.p.A. is glad to inform you that the Automated Endoscope Disinfectors models **EW 1** and **EW 1 S** - used in combination with **SteelcoXide chemicals** - can claim a **double conformity** to:

- High Level Disinfection (HLD) according to EN 15883-4
- Liquid Sterilization (LS) according to EN 14937

within the same unit, only selecting the cycle desired in the display of the machine.

This is a step ahead in enhancing assurance in endoscopes reprocessing giving the possibility to run the right cycle according with the risk analysis assessed for the specific procedure, endoscope and patient history.

#### Hydrogen peroxide advantages compared to peracetic acid:

- Longer shelf life improves products stock and warehouse management
- Less odour
- Excellent material compatibilityz



Steelco Chemicals



EW 1



EW 1 S



Follow us

